

## Adding a Differential Thermistor Sensor Model 023-0319 in Alliance


This update contain some missing diagrams and a couple new sensor model but it is mainly intended to add the 023-0319 sensor model to monitor any CPC sensor connected to a CPC AI board model 810-3010 or 810-3061

You can't create a user defined model with the configure sensor model window for this model. Use the DatabaseEditor or replace the Sensor Models Table with a newer one.

To update the site:

1. Transfer the file "Sensor Models Update for MTA51.zip"
2. Close MT Alliance
3. Make a Backup
4. Open the zip file and extract the bmp files in the "C:\Alliance\Images\Sensor Models" folder
5. Then extract the Sensor Models.\* files in the "C:\Alliance\Shared data\Config" folder
6. Make a Backup a second time now or after you make change before you leave the site
7. Start MT Alliance

To monitor a CPC 10k thermistor used by a CPC AI board you can drop this type of sensor (or change an existing TEMPERATURE sensor) to the CPC 023-0319 model



Diff Input diagram Type2.bmp

